

## Z 14 Video / Digital Zoom Inspection System / Moticam 480P Camera

With a 14:1 zoom ratio, motorised zoom lens & interchangeable objective, the Z 14's full colour high resolution image zoom range provides 3.4x – 48x (0.5x) – 6.9x – 96x (1.0x) & 10x – 144x (1.5x) when displayed on a 19" video monitor. Full colour images may be captured to a PC with USB 2 port.

Z 14 Bench System



Z 14 Video Zoom System Options

- 14 – 1 Motorised Zoom
- Handset Control
- S – VHS & Phone Output
- USB 2 Digital Image Capture
- 0.5x Objective (WD = 172mm)
- 1.0x Objective (WD = 75mm)
- 1.5x Objective (WD = 40mm)
- Bench Stands
- Articulated Stands
- Long Arm Industrial Stands
- Scanning Stages
- Locking Stages
- Gimble Rotating Stage (SMD)
- Please Contact GDC For Further Details

Z 14 Scanning Stage Example



Moticam 480P



- Live Video / Digital USB Output
- Manual White Balance
- Images 2.0 Software
- Camera Image Capture To PC
- Measurement & Annotation
- Adapts To All Microscopes
- Options Include:-
- Macro Zoom Lenses
- Copy Stands

- **Purchase the Z 14 Video Zoom with confidence from GDC.**
- **Moticam 480P will fit any microscope – Live Colour Video & Digital Image Capture to a PC**
- **GDC import & distribute Motic Microscopy products – fast & friendly support.**
- **GDC also supply macro stands, macro lenses & gooseneck arm stands for wide field applications.**
- **GDC provide local support & if required on site product / software installation & training.**