

Moticam 2.0 – 5.0 Megapixel Microscope Digital Camera Kits

Images 2.0 control software is included with every camera kit, user friendly & comprehensive Images 2.0 requires a Windows XP or Vista PC platform with a USB 2 port. Images 2.0 software features include, Camera control & Image Capture, Measurement & Annotation, Scale bar & Cross line generation. A Calibration slide, eyepiece adapters & macro lens are also included as standard. Upgrade Option - Advanced 3.2 Image Analysis / EDF Software

Moticam 2000



2.0 Megapixel

- 1600 x 1200 Pixels
- USB 2 / Twain Compatible
- C Mount Camera
- Eyepiece Mounts Included
- Images 2.0 Software
- Calibration Slide
- Contact GDC For Current List Price


3.0 Megapixel

- 2048 x 1536 Pixels
- USB 2 / Twain Compatible
- C Mount Camera
- Images 2.0 Software
- Eyepiece Mounting Kit
- Calibration Slide
- Contact GDC For Current List Price

Moticam 2300



Moticam 2500



5.0 Megapixel

- 2048 x 1536 Pixels
- USB 2 / Twain Compatible
- C Mount Camera
- Images 2.0 Software
- Eyepiece Mounting Kit
- Calibration Slide
- Contact GDC For Current List Price

- Buy Moticam Digital kits with confidence from GDC.
- GDC import & distribute Motic Microscopy products.
- GDC also supply macro lenses & Copy stands, gooseneck arms for macro imaging applications.
- GDC provide local support & if required on site software installation & training.