

Moticam Digital Camera Kits

Motic Images 2.0 control software CD is included within every camera kit, user friendly & comprehensive Images 2.0 requires a Windows XP or Vista PC platform with a USB 2 port. Software features include - Camera control & Image Capture, Various Measurement & Annotation functions, Scale bar & Cross line generation. A Calibration slide, eyepiece adapters & macro lens are also included with every kit as standard.

Moticam 480P



Live Video & Digital Capture Output

- Live Video Output – S-VHS / Phono output
- USB 2 digital Image Capture
- Manual White Balance
- Excellent Colour Reproduction
- Images 2.0 Software
- Calibration Slide
- Eyepiece mounts / AC Adapter & Leads

Low Cost Digital Camera Kit For Microscopes

- Entry Level Digital Microscopy Camera Kit
- 1280 x 1024 Pixels
- Images 2.0 Software
- Eyepiece Mounting Kit
- Calibration Slide
- USB 2 Cable

Moticam 1000



Moticam 353 / 483



Gooseneck Mounted Camera

- Gooseneck camera with fixed macro lens
- Live Video & Digital USB 2 Output (483)
- Digital Only 800 x 600 pixels (353)
- Microscope Eyepiece Adapters
- Optional CCTV Macro Zoom Lens
- Images 2.0 Software
- AC Adapter & Leads
- Carry Case

- **Buy Moticam Digital kits with confidence from GDC – fast & friendly support.**
- **GDC Microscopes import & distribute Motic Microscopy products.**
- **Moticam Digital camera kits installed UK user base runs into the thousands.**
- **Buy with confidence - GDC Microscopes provide local support & on site installation if required.**