

MOTIC EDUCATION MICROSCOPES – 18 / 28 SERIES SELECTION

Motic 2801 - LED



- Monocular head, inclined 45° and 360° rotating
- Widefield eyepiece WF10X/18mm
- Quadruple nosepiece
- Achromatic objective EA 4X, 10X, 40X S
- Coaxial coarse and fine focusing adjustment
- Rack & pinion focusable 1.25 N.A. Abbe condenser
- Iris diaphragm with filter holder
- LED illumination 20mA, 3.5V, 70mW with intensity control
- Main supply 220V-240V, charger included (CE)

Motic 2820-LED

Binocular head, inclined 45° and 360° rotating
 Widefield eyepieces WF10X/18mm
 Quadruple nosepiece
 Achromatic objective EA 4X, 10X, 40X S
 Coaxial coarse and fine focusing adjustment
 Built in low position coaxial mechanical stage
 Rack & pinion focusable 1.25 N.A. Abbe condenser
 Iris diaphragm with filter holder
 LED illumination 20mA, 3.5V, 70mW with intensity control
 Main supply 220V-240V, charger included (CE)



PM - 1805



- Monocular head, inclined 45° and 360° rotating
- Widefield eyepieces WF10X/18mm (plain cross hair)
- Quadruple nosepiece
- Achromatic objective EA 4X, 10X, 40X S
- Separate coarse and fine focusing adjustment
- Working table rotating from 0° - 360° with graduation to 1°
- Focusable 1.25 N.A. Abbe condenser
- Iris diaphragm with filter holder
- Halogen illumination 12V/20W with intensity control
- Main supply 220V-240V (CE)
- With simple analyser on type and swing-in Polariser

Summary

- **Purchase Motic education microscopes with confidence from GDC.**
- **GDC import & distribute Motic Microscopy products – competitive price guaranteed.**
- **GDC provide local support & if required on site installation & training.**