

MOTIC EDUCATION MICROSCOPES – STEREO

ST-39C-N9GO



20x & 40x Magnification

- 45° Binocular head, 360° rotating
- Widefield eyepieces WF10X/20mm
- Objective turret magnification changer (2X, 4X)
- Incident halogen light 12V/10W
- Transmitted tungsten light 12V/10W
- Main supply 220V-240V (CE)

10x & 30x Magnification

- 45° Binocular head, 360° rotating
- Widefield eyepieces WF10X/20mm
- 3:1 zoom ratio; magnification range: 1X-3X
- Incident halogen light 12V/10W
- Transmitted halogen light 12V/10W
- Main supply 220V-240V (CE)

ST-39Z-N9GO



DM-39C-N9GO-A



Motic Digital Stereo Microscope

- Widefield eyepiece WF10X/20mm on binocular stereo head
- Built-in digital imaging chip
- 1280x1024 pixels with USB2.0 output
- Turret-type objectives 2X, 4X
- Fixed-arm type stand
- Black, white and frosted glass stage plates
- Incident halogen illumination 12V/10W
- Transmitted tungsten illumination 12V/10W
- Images 2.0 Software

Buy With Confidence From GDC Microscopes – Why?

- **Buy Motic Education microscopes with confidence from GDC.**
- **GDC import & distribute Motic Microscopy products.**
- **GDC provide local support & if required on site installation & training.**