

## Moticam Digital Camera Kits

Motic Images 2.0 control software CD is included within every camera kit, user friendly & comprehensive Images 2.0 requires a Windows XP or Vista PC platform with a USB 2 port. Software features include - Camera control & Image Capture, Various Measurement & Annotation functions, Scale bar & Cross line generation. A Calibration slide, eyepiece adapters & macro lens are also included as standard.

### Moticam 480P



#### Live Video & Digital Capture Output

- Live Video Output – S-VHS / Phono
- USB 2 digital Image Capture
- Manual White Balance
- Excellent Colour Reproduction
- Images 2.0 Software
- Calibration Slide

#### Low Cost Digital Camera Kit For Microscopes

- Entry Level Digital Microscopy Camera Kit
- 1280 x 1024 Pixels
- Images 2.0 Software
- USB 2 Connection Cable
- Eyepiece Mounting Kit
- Calibration Slide

### Moticam 1000



### Moticam 353 / 483



#### Gooseneck Mounted Camera

- Gooseneck camera with fixed macro lens
- Live Video & Digital Output (483)
- Digital Only 800 x 600 pixels (353)
- Microscope Eyepiece Adapters
- Optional CCTV Macro Zoom Lens
- Images 2.0 Software
- AC Adapter & Leads
- Carry Case

### Buy With Confidence From GDC Microscopes – Why?

- Buy Moticam Digital kits with confidence from GDC.
- GDC import & distribute Motic Microscopy products.
- GDC provide local support & if required on site installation.