

## Motic SMZ & K Series Long Arm Stands

### SMZ-168 SERIES CHOICE OF BASE STANDS

Universal Stand



SW01.09K1

- Diameter of base 380mm
- Diameter of vertical pole 25mm
- Height of pole 200mm
- Diameter of pole for focusing mount 32mm

Special Universal Stand



SW01.09K4

- Length of base 250mm
- Height of pole 373mm
- Width of base 250mm
- Diameter of pole for focusing mount 32mm

Ball Bearing Boom Stand



SP99.0004

- Diameter of vertical pole 36mm
- Diameter of pole for focusing mount 32mm

Ball Bearing Boom Stand



SP99.0006

- Diameter of vertical pole 36mm
- Diameter of pole for focusing mount 32mm

Articulating Arm Boom Stand



SP99.0003

- Diameter of vertical pole 36mm
- Diameter of pole for focusing mount 32mm

Articulating Arm Boom Stand



SP99.0005

- Diameter of vertical pole 36mm
- Diameter of pole for focusing mount 32mm

- Motic Head Holder / Focusing Arm required for above stands
- Motic Microscopes for Electronics Inspection & Reworking
- Motic Microscopes for OEM Mounting to Wire Bonders & Pull Test Machines
- Call GDC Microscopes to discuss your application & provide latest price & delivery information