

Motic BA 300 POLARISING MICROSCOPE

BA 300 POL



Standard BA 300 POL Specification

- Trinocular head Siedentopf type 30° inclined (Light split 0:100)
- Widefield high eyepoint eyepiece WFPL10X/20mm with diopter
- adjustment on both eyepieces, cross hair on one eyepiece
- Intermediate tube with rotatable analyzer 360° and Bertrand lens
- Reversed quadruple nosepiece, single positions centerable
- CCIS EF Plan Strain-free objectives: 4X, 10X, 40X(S), 60X(S)
- 360° Circular Rotating Stage
- Achromat swing-out condenser N.A. 0.90/0.13 (strain-free)
- Rotatable polariser, fixed on condenser carrier
- Blue filter, power cord, allen hexagonal key, vinyl dust cover
- Koehler illumination quartz halogen 6V/30W with intensity control
- Universal power supply 100-240V, CE
- **Plates & Compensators**
- 1101006700151 λ Tint Plate
- 1101006700161 $1/4\lambda$ Retardation Plate
- 1101006700171 Quartz Wedge
- SG069903 Empty Insertion Plate
- SP030014 Mechanical Stage
- GDC Guide Price – Please Contact GDC

Summary

Buy with confidence from the GDC partnership. Please contact us for further information on the Motic BA 300 POL microscope.

We import & distribute Motic Microscopy products, GDC partners will be pleased to organise an evaluation & provide professional & friendly local support. GDC customers are always offered on site installation & training if required.