

## Motic K Series CMO Range

Great value & performance back up by Motic Microscopes 5 year factory warranty. The K series microscopes are used in many diverse applications whenever the application demands quality optical performance, mechanical build quality & very competitive prices.

From Electronics to Entomology & Forensics to Field Studies - Choose Motic from GDC.

### DSK 500 Discussion Microscope



#### Specification

- Common main objective (CMO) Infinity optical system
- 45° Binocular head (primary and secondary units)
- Widefield eyepieces WF10X/21mm (primary and secondary units)
- Dual observation optical bridge with pointer on primary unit
- Standard 1X Achromat objective
- 5-position step magnification changer
- Magnification stops: 6.4X, 10X, 16X, 25X and 40X
- 89mm working distance
- Specialised incident large working area stand

### K Series



#### Bench Stereo Microscope with Incident & Transmitted Illumination

- Common Main objective high resolution design with 1x Plan HR lens
- 110mm Working Distance
- Binocular stereo with 6x, 12x, 25x & 50x Magnification
- Bench stand with halogen or LED incident & transmitted illumination
- Various Long arm Industrial inspection stands
- Fibre optic & LED illumination options
- Camera Port Adapter
- K series design identical to the Wild M5 stereo - components are interchangeable.
- K Series Bench Microscope with integral Illumination (Pictured Left)
- Please contact GDC for current price guide

## Summary

- **Buy with confidence from GDC,**
- **Motic Europe appointed importer & distributor of Motic Microscopy products.**
- **GDC Microscopes provide local customer support we will be happy to assist you with on site installation & product training if required.**