

Motic PSM 1000 Probe Station Optics

PSM 1000



Standard PSM 1000 Specification

- Trinocular head (Optical split 50:50/0:100)
- 3-lens changeover turret (1XUV, 1XIR, 2XVIS)
- Probe Station Focusing Block-20kgs rated
- Adjustable/centerable forward facing nosepiece
- One fixed and three adjustable positions
- Laser capabilities (532nm-1064nm)
- Laser safety kit (shims and laser safety pin)
- SG01S102C3 Plan Apo objective ELWD 2X N.A. 0.055, WD 34mm
- SG01S162C3 Plan Apo objective ELWD 5X N.A. 0.055, WD 34mm
- SG01S022C3 Plan Apo objective ELWD 10X N.A. 0.280, WD 33.5mm
- SG01S032C3 Plan Apo objective ELWD 20X N.A. 0.420, WD 20mm
- SG01S122C3 Plan Apo objective ELWD 50X N.A. 0.550, WD 13mm
- SG01S062C3A Plan Apo objective ELWD 100X High Numerical Aperture N.A. 0.800, WD 3mm
- PSM 100 Guide Price – Please Contact GDC

Summary

Buy with confidence from the GDC partnership. Please contact us for further information on the PSM 1000 microscope.

We import & distribute Motic Microscopy products, we will be pleased to organise an evaluation & provide professional & friendly local support. GDC customers are also offered on site installation & training when required.