

Motic Zoom Stereo Range

Great value & performance back up by Motic Microscopes 5 year factory warranty. The SMZ series & K series microscopes are used in many diverse applications whenever the application demands quality optical performance, mechanical build quality & very competitive prices. From Electronics to Entomology & Forensics to Field Studies - Choose Motic from GDC.

SMZ 140 Series



- 10x – 40x Zoom Stereo with standard 10x eyepieces
- Binocular & Trinocular port options
- Long arm industrial stands
- OEM focus arms
- Mechanical stages
- Polarising kits
- LED & Fibre Optic Illumination

SMZ 140 Series Price Guide From £

The SMZ 168 –TL (Pictured Right) comprises a 7.5x – 50x Zoom (with 10x Eyepieces supplied as standard) with click stop zoom & Trinocular port for video / digital image capture. The large work area stand features incident & transmitted illumination with variable intensity control. Options include X / Y Mechanical Stage, Polarising kit, LED & Fibre Optic Illumination & of course Motic Digital camera kits. Long arm industrial stands also OEM mounting focus arms are also available to adapt to equipment such as wire bonders & pull test machines.

SMZ 168 Series



GDC Price Guide – SMZ 168-TL Complete £899 + VAT

K Series



Common Main objective design with 1x Plan HR lens Binocular stereo with 6x, 12x, 25x & 50x Magnifications. Camera port option & variety of accessories Are available. Bench stands that incorporate halogen or LED incident & transmitted illumination.

The design is almost identical to the Wild M5 stereo so components are interchangeable.

GDC System price guide -

Summary

Buy with confidence from GDC, we import & distribute Motic Microscopy products, we provide local support & when required on site installation.....