


INVERTED MICROSCOPES – AE20 / AE21 / AE30 / AE31

AE 20 / AE 21 INVERTED SERIES




- o Example Specification AE 21 Trinocular:-
- o Trinocular head Siedentopf type 45° inclined (light split 20:80)
- o Widefield high eyepoint eyepiece WF10X/20mm
- o Side facing quadruple nosepiece
- o CCIS E Plan Achromat objective EPL4X, LWD EPL40X
- o CCIS E Plan Achromat phase objective EPLPH10X, LWD EPLPH20X
- o ELWD N.A. 0.30 condenser with phase slider, centrable PH1, brightfield / Centering telescope
- o Plain stage
- o Halogen illumination 6V/30W with intensity control
- o 45mm Ø blue, green interference and neutral density filters
- o Hexagonal key, vinyl dust cover
- o Universal power supply 100-240V, CE


Petri Dish and Haemacytometer Holders




Ø 35mm Petri Dish Holder



Ø 54mm Petri Dish Holder



Ø 65mm Petri Dish Holder



Haemacytometer Holder

Motic AE 30 / AE31 INVERTED SERIES



- o Example Specification – AE31 Trinocular
- o Trinocular head Siedentopf type 45° inclined (light split 20:80)
- o Widefield Eyepieces WFPL10X/22mm with diopter adjustment and rubber eyecup
- o Side facing quintuple nosepiece
- o CCIS Plan Achromat objective PL4X
- o CCIS Plan Achro Phase objective PL PH10X, LWD PL PH20X, LWD PL PH40X
- o ELWD N.A. 0.30 condenser
- o Phase slider, centrable PH1, brightfield, PH3, Center Telescope
- o Plain stage, glass and metal stage inserts
- o Koehler illumination quartz halogen 6V/30W
- o 45mm Ø blue, green interference and neutral density filters
- o Hexagonal keys, vinyl dust cover
- o Universal power supply 100-240V, CE

- **Buy Motic Laboratory Microscopes With Confidence From GDC Microscopes.**
- **Great Value, Great Performance - Back Up By The 5 Year Motic Factory Warranty.**
- **GDC are Motic Europe Appointed Importers & Distributors of Motic Microscopy Products.**
- **GDC Sales Are Locally Supported & If Required GDC Provide On Site Installation & Product Training.**