

## Moticam Digital Camera Kits – Laboratory

Every digital camera kit & digital microscope includes everything necessary to get started. Images 2.0 software is provided as standard, providing user friendly software to control & capture high resolution colour images & video clips. Measurement & annotation features are also included. The Advanced 3.2 software option adds image analysis, extra depth focusing & large image capture.

### Moticam Digital Laboratory Products



### Moticam Digital Kits

USB 2 Camera Kits with Images 2.0  
Guide Price - See website for up to date prices

- Moticam 1000 1.3 Megapixel
- Moticam 2000 2.0 Megapixel
- Moticam 2300 3.0 Megapixel
- Moticam 2500 5.0 Megapixel

IEEE1394 CCD Camera Kits  
With Advanced 3.2 Software

- Moticam 3000 3.3 Megapixel
- Moticam 3000 Cooled
- Moticam 5000 5.0 Megapixel
- Moticam 5000 Cooled

### DMBA 300 Digital Laboratory Microscope

- Complete Microscope System
- 3.0 Megapixel Digital Head
- Images 2.0 Software
- USB 2 Output
- 4x, 10x, 40x, 100x Objectives
- Low position Coaxial Stage
- Focusable & Centerable Condenser
- 6v / 30watt Quartz Halogen Koehler
- Guide Price – See website for up to date prices

### DMBA 300 Digital Microscope



- **Buy Moticam Digital Laboratory Products with confidence from GDC.**
- **GDC import & distribute Motic Digital Microscopy products.**
- **GDC also supply macro stands, gooseneck arms for macro imaging applications.**
- **GDC provide local support & if required on site product / software installation & training.**