

## Laboratory Microscopes – BA 400 Fluorescence / Multi Head



- Example BA 400 + Epi Fluorescence Specification:-
- **EF-UPR Epi-Fluorescence attachment**
- AZ1000022A Complete Epi-Fluorescence Attachment with
- AZ1000023 Epi-Fluorescence attachment with filter cassette
- SE060110 Mercury lamp HG-100W
- AZ1000040 Lamphouse HBO 100W
- AZ1000041 Collector lens
- SE020035 Mercury lamp socket for HBO 100W
- SP120066 Starter Unit HBO 100W
- AZ1000042 Lamp centering tool
- SE030021 VDE power cord
- SE030064 UK power cord
- Various Fluorescence Filter Sets Available



### Motic BA 400 DISCUSSION MULTI HEAD

- Example Specification
- **BAT-400-MVH5 Multi Observation System**
- BA-400 Trinocular Standard System (primary unit)(light split 20/80)
- Koehler illumination quartz halogen 6V/30W with external lamp house
- 4x Binocular head Siedentopf type 30°
- 10x Widefield high eyepoint eyepieces
- Multi observation optical bridge (2 ) supporting 2 secondary units
- LED pointer on primary unit
- Adjustable height support apparatus
- Universal power supply 100-240V, CE
- 2 Head & 3 Head Versions Available

- **Buy Motic Laboratory Microscopes With Confidence From GDC Microscopes.**
- **Great Value, Great Performance - Back Up By A 5 Year Motic Factory Warranty.**
- **GDC are Motic Europe Appointed Importers & Distributors of Motic Microscopy Products.**
- **GDC Sales Are Locally Supported & If Required GDC Provide On Site Installation & Product Training.**